

Face Shields

Instructions for use:

1. Do not touch the transparent sheet while wearing it.
2. The transparent sheets can be replaced as per the need.
3. Replace the Face Shield as soon as there is a splatter of bodily fluids on the sheet & dispose of the Face Shield properly.
4. Replace transparent sheets every day. The headband can be reused after proper sterilization.
5. Disinfection needs to be carried out every day. Follow the disinfection protocol

Disinfection protocol for End user: Soap/ Alcohol based sanitizers

Your face shield has 2 components

- a) A Reusable headband with elastic
- (b) A transparent sheet attached to the reusable headband.

While Using soap

1. Take a clean cloth and dip it in a soapy water.
2. Rub the headband (both sides along with the elastic) thoroughly with the cloth and allow the soap to remain on the surface for 30 seconds.
3. Wipe the soap foam from the headband with another clean cloth
4. Allow to dry before using again

While Alcohol based cleaning

1. Use a 70% isopropyl alcohol wipe to rub the surface of the headband on both the sides and in gaps and crevices. Also rub along the elastic band with the same or a new alcohol wipe.
2. In case wipes are not available use a tissue paper. Put a coin size drop of gel-based hand sanitizers (61-70% alcohol content) on the top surface of headband. Wipe it with a clean tissue paper. Repeat the same procedure for the other side of the headband and the elastic.

In case of problems:

Call – Sayali: +91-8788554002